


Neuromonitoring Analysis (NMA®) Faster reliable diagnosis.



- Innovative TCD analysis
- Screening tool for cardiovascular physiology
- Especially for the use in ICU and OR
- Advanced TCD parameter
- Faster reliable interpretation of the TCD signal
- Easy differentiation of findings
- Only available for DWL® systems




DWL Neuromonitoring Analysis (NMA[®])


The new TCD analysis software for cardiovascular physiology enables the differentiation of pathological and non-pathological findings with the aid of the Neuromon B.V. Doppler parameters and can be used, among other things, in surgery, anaesthesia and intensive care.

The use of this innovative screening software allows physicians and healthcare professionals to quickly and reliably interpret the TCD signal in complex clinical situations, providing valuable information for further disease assessment and therapy control.

Common Indices




New Parameters (Indices)




Acceleration
(maximal rate of change)
in FV per sec

Example case: pre- and postoperative examination – patient 57 years with 70% ICA stenosis




Before surgery



After surgery

Example case: therapy control medical treatment – patient with sepsis.



The company Neuromon B.V. has further developed the existing TCD indices for advanced TCD parameterization, whereupon appropriate software for the NMA[®] module has been specially designed for DWL[®] devices.

These “new parameters” are based on the theory of arterial acceleration, which states that the arterial vasculature is not just a passive conduit, but brings the energy into the pressure wave of the heart. This pressure wave is expanded within the smooth muscle cells of the arterial walls and spreads as a peristaltic wave along the branches of the Arterial Tree.

The available heart beats are combined in 10s intervals and the following parameters are determined:

- acc** acceleration, the maximal change in flow velocity at stroke onset.
- sys1** the maximal flow velocity reached during early systole (first systolic peak)
- sys2** the maximal flow velocity reached during late systole (second systolic peak)
- dias@560** the diastolic flow velocity at fixed time after stroke onset (560 ms)
- HR** mean heart rate over the 10s interval

Parameter	Origin	1st order dependence	2nd order dependence	3rd order dependence
acc (sys1)	phasic myogenic response	smooth muscle contractility		
		aorta pressure		
sys2	ventricular ejection	stroke volume	diastolic filling time	heart rate
		blood distribution	ventricular ejection fraction	heart contractility
			cerebro-vascular resistance	metabolic activity
			peripheral vascular resistance	sympathetic vasomotor tone metabolic activity
dias@560	aorta pressure	cardiac output	preload to the heart	venous capacity total blood volume
		aorta stiffness		
		total vascular resistance		
HR	cardiac innervation	aorta and carotid baro-receptors	dynamic aorta pressure	heart contractility aorta pressure
		atrial baro-receptors	venous pressure	total vascular resistance
				venous capacity
				total blood volume cardiac output



Authorized DWL® Distributor

DWL USA, Inc.
4062 Flying C Rd.
Suite 37
Cameron Park
CA 95682
USA

Phone 760-753-7319
sales@dwl.us
www.dwl.us

Neuro-Monitoring-Analysis (NMA®)

Referenzen / References:

Schaafsma, A. (2012). Improved Parameterization of the Transcranial Doppler Signal. *Ultrasound in Medicine and Biology*, 38(8), 1451-1459.

Schaafsma A. Harvey with a modern twist: How and why conducting arteries amplify the pressure wave originating from the heart. *Med Hypotheses*. 2014;82(5):589-594.

De Goede, A. A., Loef, B. G. Reidinga, A. C., & Schaafsma, A. (2017).

Fluid Resuscitation in Septic Patients Improves Systolic but Not Diastolic Middle Cerebral Artery Flow Velocity, 43(11), 2591-2600.

Manufacturer: Compumedics Germany GmbH,
Josef-Schuetzler-Strasse 2, 78224 Singen / Germany, www.dwl.de
Compumedics Germany GmbH is certified in accordance with EN ISO 13485.
Multi-Dop, EZ-Dop, Doppler-Box, DiaMon, NMA, DWL and the DWL logo are

trademarks of Compumedics Germany GmbH, Compumedics and the Compumedics logo are trademarks of Compumedics Limited Australia. Windows is a trademark of Microsoft Corporation.
GTIN: QL-Upgrade NMA 04260493312462. Subject to change without notice.

In collaboration with:



www.neuromon.eu