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Manufacturer: Compumedics Germany GmbH, Josef-Schuetzler-Strasse 2, 78224 Singen, Germany

Product(s): All **re-usable, invasive DWL® ultrasound probes** delivered by Compumedics Germany GmbH and its distributors:

REF	DI/GTIN	Description	Part of set
3652	04260493310451	16MHz probe Ø 1 mm	REF 6586 (DI/GTIN 04260493310819) 16MHz probe Ø 1mm with adapter cable for patient isolation
3654	04260493310468	16MHz probe Ø 2 mm	REF 6587 (DI/GTIN 04260493310833) 16MHz probe Ø 2mm with adapter cable for patient isolation
3655	04260493310475	16MHz probe Ø 3 mm	REF 6588 (DI/GTIN 04260493310840) 16MHz probe Ø 3mm with adapter cable for patient isolation
3653	04260493310482	16MHz probe Ø 1.8 mm	REF 6589 (DI/GTIN 04260493310826) 16MHz probe Ø 1.8mm with adapter cable for patient isolation

Product labelling: REF XXXX (lasered on connector)
Serial number: SN XXXX (lasered on connector)
CE mark: CE 0123 (lasered on connector)

Traceability: The combination of product labelling and serial number forms the unique identification for traceability.


Safety class: IP 67
Classification acc. MDD: IIa

CAUTIONS	<ul style="list-style-type: none"> - In case of infectious risk patients (Creutzfeldt-Jakob), dispose of the probe adequately after use and DO NOT reuse it. - Probe is a sensitive electromechanical high-precision instrument and must be handled with the necessary care. - The reprocessing of sterilisable medical products should only be carried out by skilled personnel. - The probe is NOT designed for temperatures above 70°C (158°F). - Do NOT clean with ultrasound. - Do NOT autoclave. - Adapter cable for patient isolation (REF 6000, 6001, 6002 and 6003) is NOT sterilisable. It can be disinfected like all other accessories (see instructions for use of device).
-----------------	--

Advice and restrictions of treatment	<p>The washer-disinfectors for reprocessing must comply with the national applicable regulations and standards (e.g. DGHM-or FDA-approval or CE-marking, EN ISO 15883).</p> <p>To avoid cross-contamination, the assembly and disassembly should only be carried out by trained personnel wearing clean gloves.</p> <p>In case of visual and/or functional damage, the probe must not be used and should be sent back to the manufacturer for repair or replacement. Please ensure that proper decontamination (cleaned and sterilized) of the probe was performed before shipment and that the probe is labelled accordingly.</p> <p>Reprocessing usually has only a limited effect on the condition of the probe. The end of product life is mainly determined by wear, tear and possible damage by use and reprocessing.</p> <p>The product life of the probe is a maximum of 25 cycles for STERRAD® sterilizers (type 100NX™) and STERICOOL® sterilizers (type A110S) if the probe is used and reprocessed with appropriate care. This also applies to the protective cap.</p> <p>Reprocess the probe as soon as possible (within one hour) after use.</p>
Identification und traceability	<p>The combination of engraved marking on connector and engraved serial number on connector is unique.</p> <p>The number of reprocessing cycles cannot be traced on the probe, so please keep separate records.</p>
Accessory	For the processing and later storage of the probe, a special sterilisation basket is necessary (REF 6474). This article is not included as standard.

REPROCESSING INSTRUCTIONS

Delivery	The probes are delivered in a non-sterile condition . They must be processed before first use. Each probe is delivered with a protective cap (REF 6542) as standard.
Point of use	<p>To avoid cross-contamination, the assembly and disassembly should only be carried out by trained personnel wearing clean gloves.</p> <p>Disconnect the probe from the adapter cable. Make sure that no liquid gets into the connector.</p> <p>Remove directly after use coarse dirt (blood, mucin ...) by wiping the probe with a lint-free damp cloth. Wrap the probe in the cloth and pull the probe through the cloth from the probe connector toward the probe tip. Repeat this procedure until no coarse dirt is visible on the probe.</p>

<p>Storage and transport</p>	<p>Close immediately the probe connector with the protective cap (REF 6542), then place the probe in the sterilisation basket and secure the probe head [see Figures I to III].</p>    <p>Figure I: probe fixed in sterilisation basket</p> <p>Figure II: Connector with protective cap (REF 6542) fixed in sterilisation basket</p> <p>Figure III: Probe head fixed in sterilisation basket</p> <p>Close the lid of the sterilisation basket for safe transport of the probe.</p>						
<p>Pre-cleaning</p>	<p>If the probe is visibly contaminated, remove dirt with a lint-free damp cloth. Further pre-cleaning is not necessary. DO NOT clean with ultrasound.</p>						
<p>Preparation for decontamination</p>	<p>Ensure that the probe connector is properly closed by the protective cap (REF 6542) and that the probe is securely fixed in the sterilisation basket.</p>						
<p>Automated cleaning</p>	<p>Make sure that the probe connector is closed with the protective cap (REF 6542).</p> <p>The following steps were carried out during the validation of the automated cleaning process using washer/disinfectors unit:</p> <ol style="list-style-type: none"> 1. 2 minutes pre-cleaning with cold tap water (10°C – 25°C [50°F - 77°F]) 2. Drain 3. 5 minutes cleaning with tap water at 55°C (131°F) and an alkaline detergent* 4. Drain 5. 3 minutes rinsing with cold, desalinated water 6. Drain 7. 2 minutes rinsing with cold, desalinated water 8. Drain and end <p>Please ensure that an alkaline cleaning detergent with a proven deactivating efficacy against pathologic prion proteins is used. A 0.5 % solution of neodisher® SeptoClean (Chemische Fabrik Dr. Weigert) was used for the validation study.</p>						
<p>Disinfection</p>	<p>Make sure that the probe connector remains closed with the protective cap (REF 6542).</p> <p>The probe is not designed for temperatures above 70°C (158°F). DO NOT carry out thermal disinfection in washer/disinfectors units! As this is a thermolabile medical product, the automated disinfection must be carried out with a washer/disinfector with the capability of chemothermal disinfection.</p> <p>During the validation of the automated disinfection process by means of washer/disinfector the following steps were carried out:</p> <ol style="list-style-type: none"> 1. 1% disinfection solution thermosept® PAA base/additive (Schülke & Mayr GmbH) in tap water at 35°C (95°F) for 5 minutes 2. pumping out 3. 3 minutes intermediate rinsing with deionized cold water (10°C - 25°C [50°F - 77°F]) 4. pumping out 5. 2 minutes final rinse with deionized cold water (10°C - 25°C [50°F - 77°F]) 6. Pump out and end 						
<p>Drying</p>	<p>The probe is not designed for temperatures above 70°C (158°F). DO NOT dry with hot air in the washer/disinfector unit.</p> <p>Manual drying can be carried out with a lint-free cloth. The cavities (connector) can be dried with sterile compressed air.</p>						
<p>Testing and maintenance</p>	<p><u>Visual inspection for damage / wear</u></p> <table border="0"> <tr> <td style="padding-right: 20px;">Probe</td> <td>The coating of the probe head must not be damaged, therefore carry out a visual inspection with a magnifying glass or under a microscope if necessary and check the probe cable for damage.</td> </tr> <tr> <td>Sterilisation basket</td> <td>Perform a visual inspection for intactness (e.g. corrosion on the sterilisation basket, wear on the flexible O-rings).</td> </tr> <tr> <td>Protective cap</td> <td>Perform a visual inspection for intactness (e.g. wear on the protective cap).</td> </tr> </table>	Probe	The coating of the probe head must not be damaged, therefore carry out a visual inspection with a magnifying glass or under a microscope if necessary and check the probe cable for damage.	Sterilisation basket	Perform a visual inspection for intactness (e.g. corrosion on the sterilisation basket, wear on the flexible O-rings).	Protective cap	Perform a visual inspection for intactness (e.g. wear on the protective cap).
Probe	The coating of the probe head must not be damaged, therefore carry out a visual inspection with a magnifying glass or under a microscope if necessary and check the probe cable for damage.						
Sterilisation basket	Perform a visual inspection for intactness (e.g. corrosion on the sterilisation basket, wear on the flexible O-rings).						
Protective cap	Perform a visual inspection for intactness (e.g. wear on the protective cap).						

Preparations and packaging	<p>For sterilisation, the probe must be stored in the sterilisation basket. Remove the protective cap (REF 6542) from the connector and leave the protective cap in the sterilisation basket. Check the position of the connector and the probe head (they must each be fixed between the flexible O-rings [see Figures I to III]).</p> <p>Double wrap the closed sterilisation basket in hydrogen peroxide pouches, which comply to EN ISO 11607 series of standards. Validation was performed with Tyvek® sterilisation pouches:</p> <table border="0"> <tr> <td>R43 inner pouch</td> <td>200mm x 400mm</td> </tr> <tr> <td>R44 outer pouch</td> <td>250mm x 500mm</td> </tr> </table>	R43 inner pouch	200mm x 400mm	R44 outer pouch	250mm x 500mm						
R43 inner pouch	200mm x 400mm										
R44 outer pouch	250mm x 500mm										
Sterilisation	<p>The probe is not designed for temperatures above 70°C (158°F). Thermal steam-sterilisation must NOT be performed.</p> <p>Sterilisation can be performed using a low temperature hydrogen peroxide sterilisation process, such as the standard cycle of Getinge AB manufactured low temperature hydrogen peroxide plasma sterilizer STERICOOL® A110S. As specified by the manufacturer, this cycle corresponds to the standard cycle of the STERRAD® 100NX™ low temperature hydrogen peroxide plasma sterilizer manufactured by ASP®.</p>										
Storage	<p>The following storage conditions must be observed:</p> <table border="0"> <tr> <td>Temperature</td> <td>-40°C to +70°C (-40°F to +158°F)</td> </tr> <tr> <td>Air pressure</td> <td>500 hPa to 1060 hPa</td> </tr> <tr> <td>Relative humidity</td> <td>10 % to 100 % (no condensation)</td> </tr> </table> <p align="center">Protect from direct sunlight.</p>	Temperature	-40°C to +70°C (-40°F to +158°F)	Air pressure	500 hPa to 1060 hPa	Relative humidity	10 % to 100 % (no condensation)				
Temperature	-40°C to +70°C (-40°F to +158°F)										
Air pressure	500 hPa to 1060 hPa										
Relative humidity	10 % to 100 % (no condensation)										
Disposal	<p>Ensure that proper decontamination (cleaning and sterilisation) of the used probe has been performed prior to disposal. If this is not done, it must be disposed of as infectious hazardous waste.</p>										
Additional information	<p>When processing several instruments in one reprocessing cycle, the maximum load capacity of the respective washer/disinfector unit and steriliser unit must not be exceeded.</p>										
Reprocessing validation study information	<p>The following equipment, materials and machines have been used to validate reprocessing:</p> <table border="0"> <tr> <td>Detergent</td> <td>neodisher® SeptoClean (Chemische Fabrik Dr. Weigert GmbH & Co. KG)</td> </tr> <tr> <td>Desinfektionsmittel</td> <td>thermosept® PAA base und additive (Schülke & Mayr GmbH)</td> </tr> <tr> <td>Washer/Disinfector Unit</td> <td>Miele Professional G 7836 CD</td> </tr> <tr> <td>Two-level Rack</td> <td>Miele E327</td> </tr> <tr> <td>Steriliser Unit</td> <td>low temperature hydrogen peroxide plasma sterilizer STERICOOL® A110S</td> </tr> </table>	Detergent	neodisher® SeptoClean (Chemische Fabrik Dr. Weigert GmbH & Co. KG)	Desinfektionsmittel	thermosept® PAA base und additive (Schülke & Mayr GmbH)	Washer/Disinfector Unit	Miele Professional G 7836 CD	Two-level Rack	Miele E327	Steriliser Unit	low temperature hydrogen peroxide plasma sterilizer STERICOOL® A110S
Detergent	neodisher® SeptoClean (Chemische Fabrik Dr. Weigert GmbH & Co. KG)										
Desinfektionsmittel	thermosept® PAA base und additive (Schülke & Mayr GmbH)										
Washer/Disinfector Unit	Miele Professional G 7836 CD										
Two-level Rack	Miele E327										
Steriliser Unit	low temperature hydrogen peroxide plasma sterilizer STERICOOL® A110S										
Notes for Service / Repair	<p>Only disinfected probes are accepted for service/repair. Please ensure that proper decontamination (cleaned and sterilized) of the probe was performed before shipment and that the probe is labelled accordingly. This must be indicated on the accompanying documents.</p>										
Manufacturer	 <p>Compumedics Germany GmbH Josef-Schuetzler-Strasse 2 78224 Singen Germany t: +49 7731 79769-0 e: info@dwl.de, www.dwl.de</p>										

The above instructions have been validated as 'SUITABLE' by the medical device manufacturer for the preparation of a medical device for re-use. The reprocessor is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired results. The state of the art and often-national laws require these processes and procedures to be validated and regularly monitored.

Tyvek is a registered trademark of the company DuPont.
 STERICOOL is a registered trademark of the company Getinge AB.
 STERRAD, ASP and NX are registered trademarks of the company Advanced Sterilization Products.
 neodisher is a registered trademark of the company Chemische Fabrik Dr. Weigert GmbH & Co. KG.
 thermosept is a registered trademark of the company Schüle & Mayr GmbH.

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STERICOOL ist ein eingetragenes Warenzeichen der Firma Getinge AB.
neodisher ist ein eingetragenes Warenzeichen der Firma Chemische Fabrik Dr. Weigert GmbH & Co. KG.
thermosept ist ein eingetragenes Warenzeichen der Firma Schüle & Mayr GmbH.