


Reliable, efficient, and user-friendly.




## DWL TCD – The gold standard screening tool for diagnostics and treatment of Sickle Cell Disease (SCD)

- **Gold standard TCD**  
TCD is the gold standard for assessing cerebral blood flow velocities and predicting the risk of stroke.
- **The use of TCD for Sickle Cell Disease**  
With TCD, healthcare professionals can detect steno-occlusive lesions in large cerebral arteries and microangiopathic lesions, which can cause ischemic stroke and silent infarcts in people with SCD.
- **DWL for Sickle Cell Disease**  
The SCD specific software functions ensure reliable and accurate findings not only for diagnostics, but also for prognosis and for guiding therapy and treatment.
- **DWL specific Sickle Cell Disease examination program**  
The SCD specific examination program automatically displays the vessel sections with the highest blood flow velocities and determines the relevant findings quickly and easily.
- **DWL specific Sickle Cell Disease report**  
The SCD specific report highlights the highest mean velocity segment for fast identification of critical conditions and quick and accurate diagnosis.







Report options

Customizable target values for up to 4 different conditions

### Standard Target Values

Normal	< 170 cm/s	Unmarked (black)
Conditional	170 to 200 cm/s	Orange
Abnormal	> 200 cm/s	Red

\*Adams, et al.

Patient ID: Name:SCD\_normal Routine\_TCD DOB: 1/1/0001  
Date of examination: 6/6/2019

Probe (MHz)	Depth (cm)	Right		Unilateral		Left	
		Mean (cm/s)	Mean (cm/s)	Mean (cm/s)	Mean (cm/s)	Mean (cm/s)	Mean (cm/s)
MCA	2	30	120	15		123	12
MCA	2	30	121	19			
MCA	2	32				121	15
MCA	2	34				128	14
MCA	2	36	0	0		136	22
MCA	2	36	131	24			
MCA	2	38	145	29		166	19
MCA	2	40	150	27		171	23
MCA	2	42	159	20		162	25
MCA	2	42	153	20			
MCA	2	42	158	19			
MCA	2	42	151	24			
MCA	2	44	145	43		127	26
MCA	2	46	137	53		128	53
MCA	2	48	141	70		125	91
MCA	2	50	187	77		103	90
MCA	2	50				117	101
MCA	2	52	205	82		103	86
BIF	2	48	108	109			
BIF	2	54				113	56
tICA	2	48	34	23			
tICA	2	48	114	112			

\* manually recalculated values are highlighted in yellow  
 Normal < 170  
 Low conditional 170 - 185  
 High conditional 185 - 200  
 Abnormal ≥ 200

Note: These values are only examples.

Date: 6/6/2019 10:45:19 AM Mean: 12015  
 Probe: 2 PW PI: 0.9753.02  
 Sample Volume: 6  
 Gain(%): 19  
 Power(SPTA): 200  
 Depth(mm): 30  
 Scale(Hz): 8047  
 Label: MCA\_R  
 Filter: SD  
 Angle Correction: 0

Date: 6/6/2019 10:52:44 AM Mean: 12912  
 Probe: 2 PW PI: 0.8333.11  
 Sample Volume: 6  
 Gain(%): 25  
 Power(SPTA): 200  
 Depth(mm): 30  
 Scale(Hz): 10206  
 Label: MCA\_L  
 Filter: SD  
 Angle Correction: 0

Date: 6/6/2019 10:45:47 AM Mean: 12119  
 Probe: 2 PW PI: 0.9425.18  
 Sample Volume: 6  
 Gain(%): 19  
 Power(SPTA): 200  
 Depth(mm): 30  
 Scale(Hz): 8047  
 Label: MCA\_R  
 Filter: SD  
 Angle Correction: 0

### DWL specific SCD report includes:

- Customizable SCD report in PDF format
- Customizable target values for up to 4 different conditions (e.g., normal/low conditional/ high conditional/abnormal)
- Highlighted target values for fast identification of critical conditions
- Automatic color marking of the target values
- Options to include only highest mean values or selected recordings for simple assessment of critical conditions
- Option to include spectrograms
- Export options via DICOM and HL7



**Authorized DWL® Distributor**

DWL USA, Inc.  
4062 Flying C Rd.  
Suite 37  
Cameron Park  
CA 95682  
USA

Phone 760-753-7319  
sales@dwl.us  
www.dwl.us

# Sickle Cell Disease

## References:

What drives transcranial Doppler velocity improvement in paediatric sickle cell anaemia: analysis from the Sickle Cell Clinical Research and Intervention Program (SCCRIP) longitudinal cohort study.  
*Jeremie H. Estep, Ze Cong, Irene Agodoa, Guolian Kang, Juan Ding, M. Beth McCarville, Jane S. Hankins, Winfred C. Wang; (2022)*

Transcranial Doppler screening in Nigerian children with sickle cell disease: A 10-year longitudinal study on the SPPIBA cohort.  
*IkeOluwa A. Lagunju, Adeyemi Labaeka, Joy N. Ibeh, Adebola E. Orimadegun, Biobele J. Brown, Olugbemi O. Sodeinde; (2021)*

Primary stroke prevention for sickle cell disease in north-east Italy: The role of ethnic issues in establishing a Transcranial Doppler screening program.  
*Colombatti R., Meneghetti G., Ermani M., Pierobon M., Sainati L.; (2009)*

Big strokes in small persons.  
*Adams R. J.; (2007)*

Prevention of a first stroke by transfusions in children with sickle cell anemia and abnormal results on transcranial Doppler ultrasonography.  
*Adams, R. J., McKie, V. C., Hsu, L., Files, B., Vichinsky, E., Pegelow, C., et al.; (1998)*

*Compumedics Germany GmbH is certified in accordance with ISO 13485 and applies according to Annex II, Section 3 of the Directive 93/42/EEC on Medical Devices. Multi-Dop, EZ-Dop, Doppler-Box, DiaMon, NMA, DWL and the DWL logo are trademarks of Compumedics Germany GmbH,*

*Compumedics and the Compumedics logo are trademarks of Compumedics Limited Australia. Windows is a trademark of Microsoft Corporation. All specifications are subject to change without notice.*