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Handgehaltene & Monitoringsonden

Handheld probes & Monitoring probes

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Handheld probes

CW Mode is used for acoustic orientation and display of Doppler signals from the region of interest. In PW mode, Doppler shift frequency and FFT information results only from the sample volume selected. You have the ability of toggling 4 MHz and 8 MHz transducers between PW and CW operation.

Available with 2 and 2.9 meter cable length.

2 MHz PW

Used in transcranial examinations and Monitoring applications to insonate the main arteries of the brain.

Insonation depth 30-150 mm.

1 MHz PW

To examine patients with insufficient temporal windows.

Insonation depth 30-150 mm.

2 + 2.5 MHz PW¹

Used in transcranial examinations for differentiation of solid and gaseous emboli insonating the main arteries of the brain. Insonation depth 30-150 mm.

4 MHz PW & CW

Used in extracranial (carotid arteries) and peripheral examinations. Insonation depth 12-30 mm.

8 MHz PW & CW

To insonate the peripheral arteries and veins. Insonation depth 6-20 mm.

5-12 MHz linear array^{1, 2, 3}

For carotis duplex in B-mode, Color-Doppler- and Triplex-Mode.

Monitoring Sonden / Monitoring probes

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Monitoring probes

Available with click&stay or screw topped options. The benefit of the click&stay fixation lies in the ability to loosen the probe by applying slight pressure on the outer fixation ring and hence to realign the probe particularly, quickly and simply. Available with 2.5 meter cable length. Optional with 3 and 5 meter extension cable (except for MRI).

2 MHz PW click&stay / screw topped

Used in transcranial examinations to insonate the main arteries of the brain. Insonation depth 30-150 mm.

1 MHz PW click&stay / screw topped

To examine patients with hyperostosis and those with complicated windows. Insonation depth 30-150 mm.

4 MHz PW click&stay / screw topped

Used in intracranial examinations. Insonation depth 12-30 mm.

2 + 2.5 MHz PW click&stay / screw topped¹

Used in transcranial examinations for differentiation of solid and gaseous emboli insonating the main arteries of the brain. Insonation depth 30-150 mm.

2 MHz MRSafe click&stay / screw topped

MRSafe probes (Magnetic Resonance Safe) do not contain any magnetic components. Available with 3 or 5 meter extension cable.

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16 MHz probes

16 MHz reusable probes

For non-invasive and invasive measurements of blood flow conditions in the microvascular and endoscopic region.

Insonation depth 0.5 mm to 5.0 mm.

Standard

- Available in 1, 2 and 3 mm diameters
- Cable length 2.5 m

Endoscopic

- Rigid cable
- Available in 1.8 mm diameter
- Cable length 2.5 m

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16 MHz probe navigator²

For safe operation of 16 MHz (preferable 1 mm) probes during surgical procedures.

Consisting of a handle and a mounted shaft made of stainless steel.

Probe navigator neuro set

- Handle
- 4 Bayonet shafts (2 long / 2 short)
- Storage tray
- Sterilization basket

Probe navigator combination set

- Handle
- 2 Bayonet shafts (long / short)
- 2 Straight shafts (long / short)
- Storage tray
- Sterilization basket

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Compumedics Germany GmbH is certified in accordance with ISO 13485 and applies according to Annex II, Section 3 of the Directive 93/42/EEC on Medical Devices. Multi-Dop, EZ-Dop, Doppler-Box, NMA, DiaMon, DWL and the DWL logo are trademarks of Compumedics Germany GmbH, Compumedics and the Compumedics logo are trademarks of Compumedics Limited Australia. Windows is a trademark of Microsoft Corporation. All specifications are subject to change without notice.

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¹ Not available in USA

² Merchandise

³ Only expandable with appropriate hardware

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Compumedics Ltd. Australia / www.compumedics.com