


www.dwl.us  
760.753.7319

**DWL**  
USA



## Doppler-Box™X1 Top performance for routine diagnostics.


- Digital DWL® Doppler technology with superior signal and image quality
- Especially for routine diagnostics
- Measurement range up to 1200 cm/s
- Multi-depth display with 9 Doppler spectra windows simultaneously
- For connecting to any external computer
- DICOM and HL7




# Doppler-Box™ X1

Dimension L x W x H	11 x 4 x 7 inch
Weight	5.5 lbs
Doppler frequencies (MHz)	PW: 1 / 2 / 16 PW & CW: 4 / 8
Doppler channels	1
Spectra windows	1-9
Doppler-M-Mode	●
Gates	max. 8000
Audio playback	●
Report generator	●
Display	Depending on PC
Resolution	Depending on PC
USB remote control	optional
Hard disk	Depending on PC
USB ports	Depending on PC
Ethernet	10 / 100 / 1000 MBit
Graphic card interfaces	Depending on PC
Minimum requirements PC	Information from distributor available
Analog outputs	4
Minimum Operating system	Windows® 8.1 Pro
Offline reader	●
Remote support	●
DWL® database	●
Archiving internal	●
Archiving external	●
Data export	UTF-8
DICOM / HL7	optional
Wireless	Depending on PC

## Handheld probes



## Microprobes



## Software

- Routine

## Basic configuration

Digital DWL® Doppler device with routine software incl. Doppler-M-Mode, probe holder

## Additional accessories

PC with operating system Windows®<sup>4</sup>, remote control, system cart<sup>4</sup>, configurable foot switch, analog output cable, system bag

## Certification



## Classification

Transcranial Doppler Device: MDD Class II a

Authorized DWL® Distributor



**Manufacturer:** Compumedics Germany GmbH, Josef-Schuetzler-Strasse 2, 78224 Singen / Germany, [www.dwl.de](http://www.dwl.de). Compumedics Germany GmbH is certified in accordance with ISO 13485 and applies according to Annex II, Section 3 of the Directive 93/42/EEC on Medical Devices. Multi-Dop, EZ-Dop, Doppler-Box, DiaMon, NMA, DWL and the DWL logo are trademarks of Compumedics Germany GmbH, Compumedics and the Compumedics logo are trademarks of Compumedics Limited Australia. Windows is a trademark of Microsoft Corporation. GTIN 04260493313285. Subject to change without notice. <sup>4</sup> Merchandise